

FISHERBRAND™ STERILE POLYSTYRENE DISPOSABLE SEROLOGICAL PIPETS WITH MAGNIFIER STRIPE - 25 ML

Fisherbrand™ Sterile Polystyrene Disposable Serological Pipets with Magnifier Stripe

The high-contrast black markings make the volumes easier to read, for improved optical clarity. Manufactured with optically clear polystyrene for accurate measurement. Fisherbrand™ Sterile Polystyrene Disposable Serological Pipets with Magnifier Stripe reduce reading errors and are designed for tissue culture work.

- Color coded magnifier stripe runs lengthwise along pipet, making meniscus obvious at arm's length – reducing reading errors
- Sterile and nonwetable pipets
- Offer a tolerance of $\pm 2\%$ at full volume and have plug-in mouth end
- 1, 2, 5, 10, 25, and 50mL pipets have ascending and descending graduations
- Improved tip design on 25 and 50mL pipets increases the effect of surface tension to minimize dripping
- Color coding on each pipet and on inner bags for multipack pipets make it easy to spot the size needed
- Each pipet lot is tested for sterility, accuracy and endotoxin level

Cat. #: Fisher Scientific 13-678-11

Certifications/Compliance: Non-pyrogenic, USP XXII, USP XXII Bacterial Endotoxins Test (85); DNase/RNase Free; BSE/TSE Free

Sterility: Sterile

Graduations: 0.2 mL

Quantity: 200/Cs., 25/Pk., 8 Packs/Cs.

Packaging Type: Individually Wrapped, paper-plastic wrap

Volume (Metric): 25 mL

Disposable: Single-use

Material: Polystyrene

Graduated: Yes

Negative Graduations: 10 mL

Tolerance: ± 0.5 mL

Color: Red

Type: Serological Pipet

Disposable: Single-use



Note: This product was the same Fisherbrand™ , but with 50 %discount