

HRP-GOAT ANTI-RABBIT IGG

HRP-Goat Rabbit IgG

货号: 29002

规格: 100 μ L

Target: Rabbit

Conjugate: Horseradish Peroxidase

描述: HRP-Goat Rabbit IgG

背景: Whole IgG antibodies are isolated as intact molecules from antisera by immunoaffinity chromatography. They have an Fc portion and two antigen binding Fab portions joined together by disulfide bonds and therefore they are divalent. The average molecular weight is reported to be about 160 kDa. The whole IgG form of antibodies is suitable for the majority of immunodetection procedures and is the most cost effective

经过测试的应用: Elisa, IP, IHC

Dilution:

1:500 – 1:5,000 for immunohisto/cytochemistry

1:5,000 – 1:100,000 for ELISA and Western blotting with chromogenic substrates

1:10,000 – 1:200,000 for Western blotting with ECL substrates

Molecular Weight: 160 kD

种属反应性: Goat

形式: Liquid

克隆性: Polyclonal

成分: 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6

Stabilizer: 15 mg/ml Bovine Serum Albumin (IgG-Free, Protease-Free)

储存条件: Store at 4 \times , 6 months; -20 \times , 2 years