

CORNING® 500 ML EASY GRIP POLYSTYRENE STORAGE BOTTLES WITH 45 MM CAPS

Corning® 500 mL Easy Grip Polystyrene Storage Bottles with 45 mm Caps

These 500 mL Corning® disposable polystyrene bottles are ideal for storage of media, buffers and other aqueous solutions.

- Easy grip sides facilitate handling
- 45 mm polyethylene plug seal cap provides an airtight seal
- Sterile, nonpyrogenic

Cat. #: Corning® 430282

Qty./Pk: 2 / Pk

Qty./Cs: 24 / Cs

Brand: Corning®

Capacity: 500 mL

Neck Style: Wide Mouth

Top Style: Threaded

Cap Style: Plug seal

Cap Material: High density polyethylene

Bottle Shape: Cylindrical

Cap Color: Orange

Bottle Feature: Graduated

Feature: Disposable

Sterile: Yes

Application: Pharmaceutical Quality Control and Life Sciences



Note: This product was Corning® Brand but 50% off Corning list price.