

## **THERMO SCIENTIFIC™ NALGENE™ RAPID-FLOW™ STERILE FILTER STORAGE BOTTLES - 1 L**

### **Thermo Scientific™ Nalgene™ Rapid-Flow™ Sterile Filter Storage Bottles**

Save time and minimize contamination with gamma-sterilized Thermo Scientific™ Nalgene™ Rapid-Flow™ Filter Storage Bottles. Use with Nalgene Bottle-Top filters or as replacement/extra receptacles for filter units. Excellent as single-use storage containers.

- Headspace allows for extra serum or other additives to be added to the filtrate
- Leakproof† screw caps to eliminate pH shift in bottle for longer life of stored media
- Sterilized by gamma radiation and individually bagged for a 5-year sterile shelf life
- Non-pyrogenic and non-cytotoxic

**Cat. #:** Thermo Scientific™ 4551000; Fisher Scientific 09-740-25F

**Material:** Virgin Clear Polystyrene Bottle, White Pigmented Virgin HDPE Closure

**Model:** Wide

**No. per Case:** 12

**Volume (Metric):** 1 L

**Sterility:** Sterile

**Mouth:** Wide

**No. per Pack:** 12

**Fit for Purpose:** Sterile

