

## **THERMO SCIENTIFIC™ NALGENE™ RAPID-FLOW™ STERILE DISPOSABLE FILTER UNITS WITH PES, CN, SFCA OR NYLON MEMBRANES - 250 ML**

### **Thermo Scientific™ Nalgene™ Rapid-Flow™ Sterile Disposable Filter Units with PES, CN, SFCA or Nylon Membranes**

Increase filtration efficiency. PES membranes and exclusive Rapid-Flow design provide fast flow rates and high throughput while defending against cell culture contamination.

- Devices feature the fastest flow rate and lowest clogging of comparable alternatives
- 0.1um PES available for mycoplasma protection
- 0.2um PES devices are stem cell and mouse embryonic assay (MEA) tested
- Non-pyrogenic and endotoxin free - Rapid-Flow filter units meet the most rigorous testing standards for any vacuum filtration unit.

**Cat. #:** Thermo 5680020

**ACT Label:** Yes

**Material (Membrane):** PES

**Sterility:** Sterile

**No. per Pack:** 12

**Volume (Metric) Receiver:** 250 mL

**Diameter (Metric) Membrane:** 50 mm

**Disposable:** Single-use

**Pore Size:** 0.2 µm

**No. per Case:** 12

**Type:** Filter Unit

**Volume (Metric) Upper:** 250 mL

**Fit for Purpose:** Sterile



**Note:** This product was the same Thermo Scientific™ brand, but with 50 %discount