

AMHR2 (DM186) RABBIT MAB

目录: 28529

产品名称: AMHR2(DM186) Rabbit Monoclonal Antibody

基因符号: AMHR; MISR2; MISRII; MRII

描述: AMHR2 antibody(DM186) Rabbit Monoclonal Antibody

背景: This gene encodes the receptor for the Mullerian hormone (AMH) which; in addition to testosterone; results in male sex differentiation. AMH and testosterone are produced in the testes by different cells and have different effects. Testosterone promotes the development of male genitalia while the binding of AMH to the encoded receptor prevents the development of the mullerian ducts into uterus and Fallopian tubes. Mutations in this gene are associated with persistent Mullerian duct syndrome type II. Alternatively spliced transcript variants encoding different isoforms have been identified.

经过测试的应用: ELISA; Flow Cyt

推荐稀释度: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human AMHR2

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with AMHR2 (DM186) on Expi293 cells transfected with human AMHR2 (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).