

## ACETYL-LYSINE RABBIT PAB

**Cat.#:** N226080

**Product Name:** Anti-Acetyl-Lysine Rabbit pAb

**Synonyms:** -

**UNIPROT ID:** -

**Background:** Acetylation of lysine residues within proteins has emerged as an important mechanism used by cells to overcome this repression.

**Immunogen:** Acetylated BSA

**Applications:** WB

**Recommended Dilutions:** WB: 1/500-1/1000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Refer to figures

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Species independent


**Conjugation:** Unconjugated

**Modification:** Acetylated

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Epigenetics and Nuclear Signaling

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Acetyl-Lysine in HeLa lysates using Pan Acetyl-Lysine antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---