

## ACETYL-SMAD2/3 (LYS19) RABBIT PAB

**Cat.#:** N225507

**Product Name:** Anti-Acetyl-Smad2/3 (Lys19) Rabbit pAb

**Synonyms:** JVI5-2; SMAD 3; hMAD-2; SMAD 2

**UNIPROT ID:** Q15796

**Background:** The protein encoded by this gene belongs to the SMAD, a family of proteins similar to the gene products of the Drosophila gene 'mothers against decapentaplegic' (Mad) and the C. elegans gene Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways.

**Immunogen:** Synthesized peptide derived from human Smad2/3 (Acetyl-Lys19)

**Applications:** IHC-P,WB

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 52 kDa; Observed MW: 60 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human,Mouse,Rat

**Conjugation:** Unconjugated

**Modification:** Acetylated

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Signal Transduction

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing