

ALDH2 RABBIT PAB

Cat.#: N225948

Product Name: Anti-ALDH2 Rabbit pAb

Synonyms: ALDH2; ALDM; Aldehyde dehydrogenase; mitochondrial; ALDH class 2; ALDH-E2; ALDH1

UNIPROT ID: P05091

Background: Alcohol metabolism; ethanol degradation; acetate from ethanol: step 2/2.

Immunogen: A synthesized peptide derived from human ALDH2

Applications: WB,IHC-P,ICC/IF,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 56 kDa; Observed MW: 56 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Tags & Cell Markers

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ALDH2 in HepG2 lysates using ALDH2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
