

ALDOA RABBIT PAB

Cat.#: N225929

Product Name: Anti-ALDOA Rabbit pAb

Synonyms: ALDOA; ALDA; Fructose-bisphosphate aldolase A; Lung cancer antigen NY-LU-1; Muscle-type aldolase

UNIPROT ID: P04075

Background: Plays a key role in glycolysis and gluconeogenesis. In addition, may also function as scaffolding protein.

Immunogen: A synthesized peptide derived from human Aldolase

Applications: WB,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 39 kDa; Observed MW: 39 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Aldolase in A549 lysates using ALDOA antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
