

ALIX RABBIT PAB

Cat.#: N225934

Product Name: Anti-ALIX Rabbit pAb

Synonyms: AIP1; ALIX; HP95; DRIP4

UNIPROT ID: Q8WUM4

Background: Class E VPS protein involved in concentration and sorting of cargo proteins of the multivesicular body (MVB) for incorporation into intraluminal vesicles (ILVs) that are generated by invagination and scission from the limiting membrane of the endosome.

Immunogen: A synthesized peptide derived from human ALIX

Applications: WB,FC

Recommended Dilutions: WB: 1/500-1/1000 FC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 96 kDa; Observed MW: 96 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ALIX in K562 lysates using ALIX antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
