

ALPHA ACTININ (1A5) MOUSE MAB

Cat.#: N261300

Product Name: Anti-alpha Actinin (1A5) Mouse Monoclonal Antibody

Synonyms: actinin 1 smooth muscle; actinin alpha 1; ACTN 1; FLJ40884

UNIPROT ID: P12814

Background: F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Immunogen: Purified recombinant protein expressed in E.coli.

Applications: IHC-P

Recommended Dilutions: IHC: 1/50-1/100

Host Species: Mouse

Clonality: Mouse Monoclonal

Clone ID: 1A5-7D1-9H10

MW: -

Isotype: IgG1

Purification: Affinity Purified

Species Reactivity: Human,Rat,Mouse


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemical analysis of paraffin-embedded Human tonsils using alpha Actinin (1A5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma Tissue using alpha Actinin (1A5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.