

ANGIOPOIETIN 2 RABBIT MAB

Cat.#: N261840

Product Name: Anti-Angiotensin 2 Rabbit Monoclonal Antibody

Synonyms: Angiotensin-2 (ANG-2)

UNIPROT ID: O15123

Background: Can induce tyrosine phosphorylation of TIE2. Binds to TIE2 receptor and counteracts blood vessel maturation/stability mediated by angiotensin-1. Its function may be context-dependent.

Immunogen: A synthetic peptide of human Angiotensin 2

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-7E6

MW: Calculated MW: 57 kDa; Observed MW: 57 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Angiotensin 2 in rat Brain lysates using Angiotensin 2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
