

## ANKRD20A1 RABBIT PAB

**Cat.#:** N225201

**Product Name:** Anti-ANKRD20A1 Rabbit pAb

**Synonyms:** ANKRD20A1; ANKRD20A; Ankyrin repeat domain-containing protein 20A1

**UNIPROT ID:** Q5TYW2

**Background:** Located in plasma membrane.

**Immunogen:** Synthesized peptide derived from the Internal region of human A20.

**Applications:** WB,ELISA

**Recommended Dilutions:** WB: 1/500-1/1000 ELISA: 1/10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 94 kDa; Observed MW: 94 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human

**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cell Biology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ANKRD20A1 in 3T3 lysates using ANKRD20A1 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---