

ANNEXIN V RABBIT MAB

Cat.#: N261683

Product Name: Anti-Annexin V Rabbit Monoclonal Antibody

Synonyms: PP4; ANX5; ENX2; RPRGL3.

UNIPROT ID: P08758

Background: This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.

Immunogen: A synthetic peptide of human Annexin V

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-3A4

MW: Calculated MW: 36 kDa; Observed MW: 36 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human, Mouse, Rat



Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of Annexin V (green) in HeLa using Annexin V antibody, and DAPI (blue) Western blot analysis of Annexin V in K562, C6, 3T3, HeLa lysates using Annexin V antibody.