

ANNEXIN VI RABBIT MAB

Cat.#: N261849

Product Name: Anti-Annexin VI Rabbit Monoclonal Antibody

Synonyms: p68; p70; ANX6; CBP68; CPB-II

UNIPROT ID: P08133

Background: May associate with CD21. May regulate the release of Ca²⁺ from intracellular stores. Miscellaneous Seems to bind one calcium ion with high affinity.

Immunogen: Recombinant protein of human Annexin-6

Applications: WB, IHC-P, IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-7G6

MW: Calculated MW: 76 kDa; Observed MW: 70 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human, Mouse, Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Annexin VI antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Annexin6 in 3T3, C6 lysates using Annexin VI antibody.