

ANNEXIN VII RABBIT MAB

Cat.#: N261848

Product Name: Anti-Annexin VII Rabbit Monoclonal Antibody

Synonyms: SNX; ANX7; SYNEXIN

UNIPROT ID: P20073

Background: Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Immunogen: Recombinant protein of human Annexin VII

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-3F6

MW: Calculated MW: 53 kDa; Observed MW: 51,47 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Annexin VII antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Annexin VII in K562, rat Brain, HeLa lysates using Annexin VII antibody.