

ARHGAP25 RABBIT MAB

Cat.#: N261869

Product Name: Anti-ARHGAP25 Rabbit Monoclonal Antibody

Synonyms: Rho GTPase activating protein 25; KAIA0053

UNIPROT ID: P42331

Background: GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

Immunogen: Recombinant protein of human ARHGAP25

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-6G6

MW: Calculated MW: 73 kDa; Observed MW: 73 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ARHGAP25 in rat brain, Raoms, Jurkat lysates using ARHGAP25 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
