

AROMATASE RABBIT MAB

Cat.#: N261870

Product Name: Anti-Aromatase Rabbit Monoclonal Antibody

Synonyms: CYP19A1; ARO1; CYAR; CYP19; Cytochrome P450 19A1; Aromatase; CYPXIX; Cytochrome P-450AROM; Estrogen synthase

UNIPROT ID: P11511

Background: Catalyzes the formation of aromatic C18 estrogens from C19 androgens.

Immunogen: A synthetic peptide of human Aromatase

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-3D8

MW: Calculated MW: 58 kDa; Observed MW: 50 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Aromatase in MCF-7, C6, 3T3, HeLa lysates using Aromatase antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
