

ATAXIN 1 RABBIT MAB

Cat.#: N261881

Product Name: Anti-Ataxin 1 Rabbit Monoclonal Antibody

Synonyms: ATX1; SCA1; D6S504E

UNIPROT ID: P54253

Background: Chromatin-binding factor that repress Notch signaling in the absence of Notch intracellular domain by acting as a CBF1 corepressor. Binds to the HEY promoter and might assist, along with NCOR2, RBPJ-mediated repression. Binds RNA in vitro. May be involved in RNA metabolism (PubMed:21475249). In concert with CIC and ATXNIL, involved in brain development .

Immunogen: Recombinant protein of human Ataxin 1

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-8H1

MW: Calculated MW: 87 kDa; Observed MW: 105 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Ataxin 1 in 3T3 lysates using Ataxin 1 antibody.