

## BBS9 RABBIT PAB

**Cat.#:** S218086

**Product Name:** Anti-BBS9 Rabbit Polyclonal Antibody

**Synonyms:** B1; D1; C18; PTHB1

**UNIPROT ID:** Q3SYG4 (Gene Accession - BC032715 )

**Background:** This gene is downregulated by parathyroid hormone in osteoblastic cells, and therefore is thought to be involved in parathyroid hormone action in bones. The exact function of this gene has not yet been determined. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jan 2017]

**Immunogen:** Fusion protein of human BBS9

**Applications:** ELISA, IHC

**Recommended Dilutions:** IHC: 150-300; ELISA: 5000-10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Isotype:** Immunogen-specific rabbit IgG


**Purification:** Antigen affinity purification

**Species Reactivity:** Human, Mouse


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

**Research Areas:** Signal Transduction, Neuroscience

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 218086(BBS9 Antibody) at a dilution of 1/160(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the fusion protein and then with 218086(Anti-BBS9 Antibody) at dilution 1/160.