

## BCL2 L14 RABBIT MAB

**Cat.#:** N263233

**Product Name:** Anti-Bcl2 L14 Rabbit Monoclonal Antibody

**Synonyms:** BCLG

**UNIPROT ID:** Q9BZR8

**Background:** Plays a role in apoptosis.

**Immunogen:** Recombinant protein of human Bcl G/BCL2L14

**Applications:** WB,ICC/IF,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R06-3D8

**MW:** Calculated MW: 37 kDa; Observed MW: 43 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cell Biology


**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of Bcl G/BCL2L14 (green) in hela using Bcl G/BCL2L14 antibody, and DAPI(blue)



Western blot analysis of Bcl G/BCL2L14 in C6, 3T3, PC-12, Raw264.7 lysates using Bcl2 L14 antibody.



Western blot analysis of Bcl G/BCL2L14 in Hela, A549, HL-60, U251, U87-MG lysates using Bcl G/BCL2L14 antibody.