

## BRPF3 RABBIT PAB

**Cat.#:** S218374

**Product Name:** Anti-BRPF3 Rabbit Polyclonal Antibody

**Synonyms:**

**UNIPROT ID:** Q9ULD4 (Gene Accession - BC117387 )

**Background:** Component of the MOZ/MORF complex which has a histone H3 acetyltransferase activity.

**Immunogen:** Fusion protein of human BRPF3

**Applications:** ELISA, IHC

**Recommended Dilutions:** IHC: 25-100; ELISA: 5000-10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Isotype:** Immunogen-specific rabbit IgG


**Purification:** Antigen affinity purification

**Species Reactivity:** Human


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

**Research Areas:** Epigenetics and Nuclear Signaling, Cell Biology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human brain tissue using 218374(BRPF3 Antibody) at a dilution of 1/40(Cytoplasm and Nucleus).



In comparison with the IHC on the left, the same paraffin-embedded Human brain tissue is first treated with the fusion protein and then with 218374(Anti-BRPF3 Antibody) at dilution 1/40.