

CACTIN RABBIT PAB

Cat.#: N225213

Product Name: Anti-Cactin Rabbit pAb

Synonyms: CACTIN C19orf29

UNIPROT ID: Q8WUQ7

Background: Involved in the regulation of innate immune response.

Immunogen: The antiserum was produced against synthesized peptide derived from the Internal region of human CACTIN. AA range:141-190

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 89 kDa; Observed MW: 95 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat

Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cactin in hepg2 lysates using Cactin antibody.