

## CATSPER4 RABBIT PAB

**Cat.#:** S217246

**Product Name:** Anti-CATSPER4 Rabbit Polyclonal Antibody

**Synonyms:**

**UNIPROT ID:** Q7RTX7 (Gene Accession - BC128139 )

**Background:** CatSper4, is a protein which in humans is encoded by the CATSPER4 gene. CatSper1 is a member of the cation channels of sperm family of protein. The four proteins in this family together form a Ca<sup>2+</sup>-permeant ion channel specific essential for the correct function of sperm cells. Voltage-gated calcium channel that plays a central role in calcium-dependent physiological responses essential for successful fertilization, such as sperm hyperactivation, acrosome reaction and chemotaxis towards the oocyte.

**Immunogen:** Fusion protein of human CATSPER4

**Applications:** ELISA, IHC

**Recommended Dilutions:** IHC: 50-200; ELISA: 1000-5000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Isotype:** Immunogen-specific rabbit IgG


**Purification:** Antigen affinity purification

**Species Reactivity:** Human


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

**Research Areas:** Neuroscience


**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing




Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 217246(CATSPER4 Antibody) at a dilution of 1/50(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the fusion protein and then with 217246(Anti-CATSPER4 Antibody) at dilution 1/50.



The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using 217246(Anti-CATSPER4 Antibody) at a dilution of 1/50.



In comparison with the IHC on the left, the same paraffin-embedded Human gastric cancer tissue is first treated with fusion protein and then with D222039(Anti-CATSPER4 Antibody) at dilution 1/50.