

## CAVEOLIN 1 RABBIT PAB

**Cat.#:** N225556

**Product Name:** Anti-Caveolin 1 Rabbit pAb

**Synonyms:** CAV1; CAV; Caveolin-1

**UNIPROT ID:** Q03135

**Background:** Caveolin-1 may act as a scaffolding protein within caveolar membranes. Interacts directly with G-protein alpha subunits and can functionally regulate their activity (By similarity). Involved in the costimulatory signal essential for T-cell receptor (TCR)-mediated T-cell activation. Its binding to DPP4 induces T-cell proliferation and NF-kappa-B activation in a T-cell receptor/CD3-dependent manner. Recruits CTNNB1 to caveolar membranes and may regulate CTNNB1-mediated signaling through the Wnt pathway.

**Immunogen:** A synthesized peptide derived from human Caveolin-1

**Applications:** WB,IHC-P,ICC/IF

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 20 kDa; Observed MW: 25 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cardiovascular


**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded mouse lung using Caveolin1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Caveolin1 in A431 lysates using Caveolin 1 antibody.



Immunofluorescence analysis of Caveolin 1 in A431 using Caveolin1 antibody.