

CD134/OX40L RECEPTER RABBIT PAB

Cat.#: N225945

Product Name: Anti-CD134/OX40L Receptor Rabbit pAb

Synonyms: TNFRSF4; TXGPIL; Tumor necrosis factor receptor superfamily member 4; ACT35 antigen; OX40L receptor; TAX transcriptionally-activated glycoprotein 1 receptor; CD antigen CD134

UNIPROT ID: P43489

Background: Receptor for TNFSF4/OX40L/GP34.

Immunogen: A synthesized peptide derived from human CD134

Applications: WB,FC

Recommended Dilutions: WB: 1/500-1/1000 FC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 29 kDa; Observed MW: 45 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD134 in Human fetal spleen lysates using CD134/OX40L Receptor antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
