

CD276 RABBIT MAB

Cat.#: N261991

Product Name: Anti-CD276 Rabbit Monoclonal Antibody

Synonyms: B7h3; B7RP-2; AU016588; 6030411F23Rik

UNIPROT ID: Q8VE98

Background: Modulates T-cell-mediated immune responses and the development of acute and chronic transplant rejection. Plays a positive regulatory role in bone formation and has a dual role in the bone-immune interface. Induces antitumor immunity as it activates both acquired and innate immunity leading to natural killer cell and CD8 T-cell dependent killing of tumor cells.

Immunogen: Recombinant protein of mouse CD276

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-9A3

MW: Calculated MW: 34 kDa; Observed MW: 90-110 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD276 in HepG2 lysates using CD276 antibody.