

CD3E(MOSUNETUZUMAB) MAB

目录: 28575

产品名称: CD3E(Mosunetuzumab) Monoclonal Antibody

基因符号: CD3e, T3E

描述: CD3E(mosunetuzumab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度 : Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CD3E

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. CD3E protein is highly expressed inside Jurkat cells. Flow cytometry analysis with 1 μ g/mL CD3E(mosunetuzumab biosimilar) mAb (28575) (Blue histogram) or isotype control mAb (Red histogram) on Jurkat cells.