

CD3G RABBIT MAB

Cat.#: N262001

Product Name: Anti-CD3G Rabbit Monoclonal Antibody

Synonyms: CD3 gamma; CD3g; T3G

UNIPROT ID: P09693

Background: The CD3 complex mediates signal transduction.

Immunogen: A synthetic peptide of human CD3G

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-4D1

MW: Calculated MW: 20 kDa; Observed MW: 20 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD3G in Jurkat lysates using CD3G antibody.



Immunocytochemistry analysis of CD3G (green) in Jurkat using CD3G antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
