

CD96 RABBIT PAB

Cat.#: N225119

Product Name: Anti-CD96 Rabbit pAb

Synonyms: CD96; T-cell surface protein tactile; Cell surface antigen CD96; T cell-activated increased late expression protein; CD96

UNIPROT ID: P40200

Background: May be involved in adhesive interactions of activated T and NK cells during the late phase of the immune response.

Immunogen: The antiserum was produced against synthesized peptide derived from the Internal region of human CD96. AA range:291-340

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 66 kDa; Observed MW: 66 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD96 in 293UV lysates using CD96 antibody.



Western blot analysis of CD96 in 293UV lysates using CD96 antibody.