

CEACAM5 (CHONGQING PRECISION) MAB

目录: 28583

产品名称: CEACAM5 (Chongqing Precision) Monoclonal Antibody

基因符号: CEACAM-5,CD66e,CEA,Meconium antigen 100

描述: CEACAM5 (Chongqing Precision) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CEACAM5

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1 μ g/mL CEACAM5 (Chongqing Precision biosimilar) mAb (28583) on Expi293 cells transfected with Human CEACAM5 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).