

CERS6 RABBIT PAB

Cat.#: S220671

Product Name: Anti-CERS6 Rabbit Polyclonal Antibody

Synonyms: CERS5; LASS6

UNIPROT ID: Q6ZMG9 (Gene Accession - NP_982288)

Background: The LASS (longevity assurance homolog) family members are highly conserved from yeasts to mammals. Six members of this family of proteins have been characterized (LASS1, LASS2, LASS3, LASS4, LASS5 and LASS6) and they all are involved in sphingolipid synthesis. LASS6 is a widely expressed 384 amino acid endoplasmic reticulum, multi-pass membrane protein. On the luminal side of the endoplasmic reticulum membrane, the N-terminal asparagine residue is glycosylated. In cells deficient for CLN9, LASS6 corrects growth and apoptosis, and increases the levels of short ceramide species, such as C14:0- and C16:0-ceramides.

Immunogen: Synthetic peptide of human CERS6

Applications: ELISA, IHC

Recommended Dilutions: IHC: 50-200; ELISA: 2000-5000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Metabolism, Cancer


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing




Immunohistochemistry analysis of paraffin embedded Human ovarian cancer tissue using 220671(CERS6 Antibody) at a dilution of 1/35(Nucleus).



In comparison with the IHC on the left, the same paraffin-embedded Human ovarian cancer tissue is first treated with the synthetic peptide and then with 220671(Anti-CERS6 Antibody) at dilution 1/35.



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using 220671(Anti-CERS6 Antibody) at a dilution of 1/35.



In comparison with the IHC on the left, the same paraffin-embedded Human brain tissue is first treated with synthetic peptide and then with D261851(Anti-CERS6 Antibody) at dilution 1/35.