

CHROMOGRANIN A RABBIT MAB

Cat.#: N262036

Product Name: Anti-Chromogranin A Rabbit Monoclonal Antibody

Synonyms: beta Granin; Chromogranin A; ChrA; SP-I; Vasostatin I; Vasostatin-2; Vasostatin II; Pancreastatin; CHGA

UNIPROT ID: P10645

Background: Pancreastatin strongly inhibits glucose induced insulin release from the pancreas.

Immunogen: A synthetic peptide of human Chromogranin A

Applications: WB,IHC-F,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-7E1

MW: Calculated MW: 51 kDa; Observed MW: 92 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Chromogranin A in PC-12 lysates using Chromogranin A antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
