

CLATHRIN HEAVY CHAIN RABBIT MAB

Cat.#: N263271

Product Name: Anti-Clathrin Heavy Chain Rabbit Monoclonal Antibody

Synonyms: Hc; CHC; CHC17; MRD56; CLH-17; CLTCL2

UNIPROT ID: Q00610

Background: Clathrin is the major protein of the polyhedral coat of coated pits and vesicles.

Immunogen: A synthetic peptide of human Clathrin heavy chain

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-6C4

MW: Calculated MW: 192 kDa; Observed MW: 192 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Clathrin heavy chain antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Clathrin heavy chain in HeLa, A549, HL-60, U251, U87-MG lysates using Clathrin Heavy Chain antibody.