

## COILIN RABBIT PAB

**Cat.#:** N225985

**Product Name:** Anti-Coilin Rabbit pAb

**Synonyms:** CLN80; COIL; Coilin; Coilin p80

**UNIPROT ID:** P38432

**Background:** Is a component of the nuclear coiled bodies (CBS) which are involved in the function or assembly/disassembly of nucleoplasmic snRNPs. During mitosis, CBS disassemble, coinciding with a mitotic-specific phosphorylation of p80 coilin.

**Immunogen:** A synthesized peptide derived from human Coilin

**Applications:** WB,FC

**Recommended Dilutions:** WB: 1/500-1/1000 FC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 63 kDa; Observed MW: 80 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Neuroscience

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Coilin in Jurkat lysates using Coilin antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---