

COLLAGEN VI ALPHA 1 RABBIT MAB

Cat.#: N262053

Product Name: Anti-Collagen VI alpha 1 Rabbit Monoclonal Antibody

Synonyms: OPLL; BTHLM1; UCHMD1

UNIPROT ID: P12109

Background: Collagen VI acts as a cell-binding protein.

Immunogen: Recombinant protein of human Collagen VI

Applications: WB,IHC-F,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-819

MW: Calculated MW: 109 kDa; Observed MW: 147 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Collagen VI in Human placenta lysates using Collagen VI alpha 1 antibody.