

COLLAGEN XIX ALPHA 1 RABBIT PAB

Cat.#: N225231

Product Name: Anti-Collagen XIX alpha 1 Rabbit pAb

Synonyms: Collagen alpha-1(XIX) chain (Collagen alpha-1(Y) chain)

UNIPROT ID: Q14993

Background: May act as a cross-bridge between fibrils and other extracellular matrix molecules.

Immunogen: Synthesized peptide derived from human COL19A1. at AA range: 161-210

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 115 kDa; Observed MW: 130 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Collagen XIX alpha 1 in HELA lysates using Collagen XIX alpha 1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
