

CRABP2 RABBIT MAB

Cat.#: N261588

Product Name: Anti-CRABP2 Rabbit Monoclonal Antibody

Synonyms: Cellular retinoic acid-binding protein II; CRABP-II

UNIPROT ID: P29373

Background: Transports retinoic acid to the nucleus.

Immunogen: Recombinant protein of human CRABP2

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-4C4

MW: Calculated MW: 16 kDa; Observed MW: 16 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of CRABP2 (green) in MCF-7 using CRABP2 antibody, and DAPI (blue)



Western blot analysis of CRABP2 in HeLa lysates using CRABP2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
