

CT A RABBIT MAB

Cat.#: N262071

Product Name: Anti-CT A Rabbit Monoclonal Antibody

Synonyms: CCT A; CCT alph; CTPCT; PCYT1; Pcyt1a

UNIPROT ID: P49585

Background: Controls phosphatidylcholine synthesis.

Immunogen: A synthetic peptide of human CT A

Applications: WB,IHC-F,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-1G4

MW: Calculated MW: 42 kDa; Observed MW: 42 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CT A in MCF-7, 3T3, HeLa lysates using CT A antibody.