

CYTOKERATIN 16 RABBIT MAB

Cat.#: N261643

Product Name: Anti-Cytokeratin 16 Rabbit Monoclonal Antibody

Synonyms: KRT16; KRT16A; Keratin; type I cytoskeletal 16; Cytokeratin-16; CK-16; Keratin-16; K16

UNIPROT ID: P08779

Background: K16 a type I cytoskeletal keratin. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. There are two types of cytoskeletal and microfibrillar keratin: type I (acidic; 40–55 kDa) [K9 to K20] and type II (neutral to basic; 56–70 kDa) [K1 to K8].

Immunogen: Recombinant protein of human Cytokeratin 16

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500–1/1000 IHC: 1/50–1/100

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-5F8

MW: Calculated MW: 51 kDa; Observed MW: 51 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at –20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Cytokeratin 16 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Cytokeratin 16 in Hela lysates using Cytokeratin 16 antibody.