

CYTOKERATIN 4 RABBIT MAB

Cat.#: N261638

Product Name: Anti-Cytokeratin 4 Rabbit Monoclonal Antibody

Synonyms: Keratin4; Cytokeratin 4; Keratin-4; CK-4; CK4; KRT14; CYK4

UNIPROT ID: P19013

Background: Cytokeratins are a subfamily of intermediate filament keratins that are characterized by a remarkable biochemical diversity, which is represented in human epithelial tissues by at least 20 different polypeptides.

Immunogen: A synthetic peptide of human Cytokeratin 4

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-2C1

MW: Calculated MW: 56 kDa; Observed MW: 56 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cytokeratin 4 in C6, HeLa lysates using Cytokeratin 4 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
