

## DAZAP1 RABBIT MAB

**Cat.#:** N262116

**Product Name:** Anti-DAZAP1 Rabbit Monoclonal Antibody

**Synonyms:** Deleted in azoospermia-associated protein 1

**UNIPROT ID:** Q96EP5

**Background:** DAZ-associated protein 1

**Immunogen:** Recombinant protein of human DAZAP1

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R05-8J7

**MW:** Calculated MW: 43 kDa; Observed MW: 43 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Epigenetics and Nuclear Signaling

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of DAZAP1 in K562, C6, 3T3, HeLa lysates using DAZAP1 antibody.