

DDX3 (6G8) MOUSE MAB

Cat.#: N261127

Product Name: Anti-DDX3 (6G8) Mouse Monoclonal Antibody

Synonyms: DBX; DDX3; HLP2; DDX14.

UNIPROT ID: O00571

Background: Multifunctional ATP-dependent RNA helicase.

Immunogen: Purified recombinant human DDX3 protein fragments expressed in E.coli.

Applications: WB,ICC/IF,IP,ChIP

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20 ChIP: 1/20

Host Species: Mouse

Clonality: Mouse Monoclonal

Clone ID: 6G8-F4-E3

MW: Calculated MW: 73 kDa; Observed MW: 73 kDa

Isotype: IgG2a

Purification: Affinity Purified

Species Reactivity: Human,Rat,Mouse,Monkey


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of DDX3 (6G8) in HeLa using DDX3 antibody.



Western blot analysis of DDX3 in HeLa, 3T3, C6, COS7, K562 and Jurkat lysates using DDX3 antibody.



Immunoprecipitation analysis of DDX3 (6G8) in HeLa lysates using DDX3 antibody.