

DDX3 MOUSE MAB

Cat.#: N261537

Product Name: Anti-DDX3 Mouse Monoclonal Antibody

Synonyms: DBX; DDX3; HLP2; DDX14.

UNIPROT ID: O00571

Background: Multifunctional ATP-dependent RNA helicase.

Immunogen: Purified recombinant human DDX3 protein fragments expressed in E.coli.

Applications: WB,IP,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

Host Species: Mouse

Clonality: Mouse Monoclonal

Clone ID: 6G8-F4-E3

MW: Calculated MW: 73 kDa; Observed MW: 73 kDa

Isotype: IgG2a

Purification: Affinity Purified

Species Reactivity: Human,Rat,Mouse,Monkey


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of DDX3 in COS7, Jurkat, C6, 3T3 and HeLa lysates using DDX3 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
