

DENN RABBIT MAB

Cat.#: N263299

Product Name: Anti-DENN Rabbit Monoclonal Antibody

Synonyms: DENN; IG20; RAB3GEP

UNIPROT ID: Q8WXG6

Background: Guanyl-nucleotide exchange factor that regulates small GTPases of the Rab family.

Immunogen: A synthetic peptide of human DENN

Applications: WB,IHC-F,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R03-8F0

MW: Calculated MW: 183 kDa; Observed MW: 220 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using DENN antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of DENN (green) in HeLa using DENN antibody, and DAPI (blue)



Western blot analysis of DENN in HeLa, A549, HL-60 lysates using DENN antibody.