

EBF3 RABBIT PAB

Cat.#: S219458

Product Name: Anti-EBF3 Rabbit Polyclonal Antibody

Synonyms: COE3; OE-2; EBF-3; HADDS; O/E-2

UNIPROT ID: Q9H4W6 (Gene Accession - BC126130)

Background: This gene encodes a member of the early B-cell factor (EBF) family of DNA binding transcription factors. EBF proteins are involved in B-cell differentiation, bone development and neurogenesis, and may also function as tumor suppressors. The encoded protein inhibits cell survival through the regulation of genes involved in cell cycle arrest and apoptosis, and aberrant methylation or deletion of this gene may play a role in multiple malignancies including glioblastoma multiforme and gastric carcinoma.

Immunogen: Fusion protein of human EBF3

Applications: ELISA, IHC

Recommended Dilutions: IHC: 100-300; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling, Neuroscience


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing




Immunohistochemistry analysis of paraffin embedded Human esophagus cancer tissue using 219458(EBF3 Antibody) at a dilution of 1/140(Nucleus).



In comparison with the IHC on the left, the same paraffin-embedded Human esophagus cancer tissue is first treated with the fusion protein and then with 219458(Anti-EBF3 Antibody) at dilution 1/140.



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 219458(Anti-EBF3 Antibody) at a dilution of 1/140.



In comparison with the IHC on the left, the same paraffin-embedded Human thyroid cancer tissue is first treated with fusion protein and then with D227142(Anti-EBF3 Antibody) at dilution 1/140.