

ERCC8 RABBIT MAB

Cat.#: N263322

Product Name: Anti-ERCC8 Rabbit Monoclonal Antibody

Synonyms: CKN1; CSA; ERCC8

UNIPROT ID: Q13216

Background: Substrate-recognition component of the CSA complex, a DCX (DDB1-CUL4-X-box) E3 ubiquitin-protein ligase complex, involved in transcription-coupled nucleotide excision repair.

Immunogen: A synthetic peptide of human ERCC8

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R03-2F3

MW: Calculated MW: 44 kDa; Observed MW: 44 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ERCC8 in A549, HL-60, U251, U87-MG lysates using ERCC8 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
