

## EYA2 RABBIT MAB

**Cat.#:** N263323

**Product Name:** Anti-EYA2 Rabbit Monoclonal Antibody

**Synonyms:** EAB1

**UNIPROT ID:** O00167

**Background:** Functions both as protein phosphatase and as transcriptional coactivator for SIX1, and probably also for SIX2, SIX4 and SIX5.

**Immunogen:** A synthetic peptide of human EYA2

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R01-6A1

**MW:** Calculated MW: 59 kDa; Observed MW: 59 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Epigenetics and Nuclear Signaling

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of EYA2 in HeLa, A549, U2OS lysates using EYA2 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---