

FACL4 RABBIT MAB

Cat.#: N263326

Product Name: Anti-FACL4 Rabbit Monoclonal Antibody

Synonyms: ACS4; FACL4; LACS4; MRX63; MRX68

UNIPROT ID: O60488

Background: Activation of long-chain fatty acids for both synthesis of cellular lipids, and degradation via beta-oxidation. Preferentially uses arachidonate and eicosapentaenoate as substrates.

Immunogen: Recombinant protein of human FACL4

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-2G7

MW: Calculated MW: 79 kDa; Observed MW: 79 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FACL4 in HeLa, A549 lysates using FACL4 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
